



## Certificate of Analysis

**Catalog #:** HFVIIa

**Lot:** Sample

**Description:** Human Factor VIIa

**Concentration:** 1.89 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.529 ml

**Activity:** 53833.00 PEU/mg

**Total protein:** 1.00 mg

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

**Storage:** < -60 C

**Expiration Date:** Sample

**PO:** Sample

**Date of Manufacture:**

The Human Factor VIIa was prepared from homogeneous Human Factor VII using Human Factor XIIa. The Factor XIIa was removed using a Corn Trypsin Inhibitor column. The FVIIa is >95% pure by SDS-PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**