



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HFVIIaGD

Lot: Sample

Description: Human Factor VIIa gla domainless

Concentration: 0.76 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.132 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 7 mM Tris-HCl/ 52.5 mM sodium chloride/ 1.4 mM EDTA/ pH 8.6/ 30% glycerol

Storage: < -60 C

Expiration Date: Sample

PO: Sample

Date of Manufacture:

The Human Factor VIIa gla domainless was prepared from purified Human Factor VIIa using Cathepsin G. The Cathepsin G is removed using affinity chromatography after digestion. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance