



Certificate of Analysis

Catalog #: HFXIIa 1212a

Lot: Sample

Description: Human Factor alpha-XIIa

Concentration: 1.13 mg/ml

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.442 ml

Activity: 72.90 PEU/mg

Total protein: 0.50 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M NaCl/ pH 5.3

Storage: < -60 C

Expiration Date: Sample

PO: Sample

Date of Manufacture:

The Human Factor XIIa was activated from homogeneous Human Factor XII using Dextran Sulfate. This Dextran Sulfate was removed after activation. Complete activation was observed on SDS-PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance