



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: HFXIIlab

Lot: SAMPLE SAMPLE

Description: Human Factor XIIa beta

Concentration: 1.17 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.085 ml

Activity: N/A

Total protein: 0.10 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M NaCl/ pH 5.3

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Human Factor XIIlab was activated from homogeneous Human Factor XII using Dextran Sulfate. This Dextran Sulfate was removed after activation. Complete activation was observed on SDS-PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance

Sample