



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HFXIa 1111a

Lot: Sample

Description: Human Factor XIa

Concentration: 1.06 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.094 ml

Activity: 193.60 U/mg

Total protein: 0.10 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M NaCl/ pH 5.3

Storage: < -60 C

Expiration Date: Sample

PO: Sample

Date of Manufacture:

The Human Factor XIa was activated from homogeneous Human Factor XI using Human Factor XIIa. This XIIa was removed using a Corn Trypsin Inhibitor column. The protein is >95% pure via SDS-PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance