



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Catalog #:** HGT

**Lot:** Sample

**Description:** Human gamma Thrombin

**Concentration:** 2.37 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.422 ml

**Activity:** N/A

**Total protein:** 1.00 mg

**Buffer:** 10 mM sodium acetate/ 0.7 M sodium chloride/ pH 6.0-HCl

**Storage:** < -60 C

**Expiration Date:** Sample

**PO:** Sample

**Date of Manufacture:**

The human Gamma-thrombin is >95% pure via SDS-PAGE and shows no reduction upon incubation with 2-mercaptoethanol. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**