



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HP-GD

Lot: Sample

Description: Human Prothrombin-gla domainless

Concentration: 8.53 mg/ml

Aliquots: 1 x 10.00 mg

Volume/aliquot: 1.172 ml

Activity: N/A

Total protein: 10.00 mg

Buffer: 50 mM Tris-HCl/ 150 mM NaCl/ pH 7.9

Storage: < -60 C

Expiration Date: Sample

PO: Sample

Date of Manufacture:

Human prothrombin with the gla domain removed. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance