



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Catalog #:** HTi

**Lot:** Sample

**Description:** Human alpha Thrombin Inactivated

**Concentration:** 3.16 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.316 ml

**Activity:** N/A

**Total protein:** 1.00 mg

**Buffer:** 50 mM sodium citrate/ 0.2 M NaCl/ 0.1% PEG-8000/ pH 6.5

**Storage:** < -60 C

**Expiration Date:** Sample

**PO:** Sample

**Date of Manufacture:**

Human alpha Thrombin treated with an irreversible inhibitor (PPACK) to block the active site. This inhibitor is removed via dialysis. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**