



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com

Certificate of Analysis

Catalog #: MFXIII

Lot: SAMPLE SAMPLE

Description: Murine Factor XIII

Concentration: 3.00 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.033 ml

Activity: 3500.00 Loewy u/mg

Total protein: 0.10 mg

Buffer: 50 mM Tris/ 100 mM NaCl/ 1 mM EDTA/ 10 u/mL aprotinin/ 20% glycerol/ pH 7.5

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The murine XIII is purified from mouse plasma and is >94% pure upon visual inspection of SDS-PAGE.

Thank You,
Quality Assurance

Sample