



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: M11a

Lot: SAMPLE SAMPLE

Description: Mouse Thrombin

Concentration: 1.28 mg/ml

Aliquots: 1 x 0.05 mg

Volume/aliquot: 0.039 ml

Activity: 2993.00 NIH u/mg

Total protein: 0.05 mg

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/ 0.1% PEG-8000/ pH 6.5

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Mouse Thrombin was activated from homogeneous Mouse Prothrombin by activation with Saw Scaled Viper Venom. The venom was removed after activation. Complete activation was observed on SDS-PAGE

Thank You,
Quality Assurance

Sample