



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human GPIIb/IIIa

Concentration: 1 mg/ml

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 08/2007

PO: N/A

Immunogen: Human Platelets

Clone ID#: 2B3A

Antibody purified by Protein G affinity chromatography from cell supernatant. Binds glycoprotein IIb/IIIa complex from platelets and megakaryocytes. Inhibits platelet aggregation and prolongs bleeding time. Binding requires calcium ion. Useful for immunofluorescence. Reconstitute to original volume with water.

Thank You,
Quality Assurance

Sample