



**1801 Commerce Drive**  
**South Bend, IN 46628**  
**Phone:(574) 288-2268**  
**Fax:(574) 288-2272**  
**Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Monoclonal antibody to Human Factor X-Light Chain

**Concentration:** 1 mg/ml (before lyophilization)

**Aliquots:** 1 x 0.50 mg

**Volume/aliquot:** 0.500 ml

**Activity:** N/A

**Total protein:** 0.50 mg

**Buffer:** 10 mM Sodium Phosphate/ 0.15 M  
NaCl/ 0.01% Sodium Azide/ pH 7.4

**Storage:** -20 °C

**Expiration Date:** 01/2010

**PO:** N/A

Immunogen: Human Factor X Light Chain

Clone ID#: 509

Antibody purified by Protein G affinity chromatography from cell supernatant. Purity determined by SDS-PAGE is >98%.

Binds human factor X and Xa in solid phase ELISA. In immunoblotting, binds factor X light chain. Bind I-125 labeled factor X. Does not bind bovine factor X. Binds factor X in presence of EDTA or Ca<sup>2+</sup>. Prolongs plasma prothrombin time (PT) and APTT. Partially inhibits prothrombin activation in prothrombinase assay.

We suggest hydrating the antibody to original concentration of 1 mg/mL.

Thank You,  
**Quality Assurance**

Sample