



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Factor XI

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 09/2007

PO: N/A

Immunogen: Human Factor XI

Clone ID#: 070

Antibody purified by Protein G affinity chromatography from cell supernatant. Binds human factor XI in solid-phase ELISA. Blots human factor XI under non reducing conditions.

We suggest hydrating the antibody to original concentration of 1 mg/mL.

Thank You,
Quality Assurance

Sample