



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Factor XII

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 02/2011

PO: N/A

Immunogen: Human Factor XII

Clone ID# GMA-140

Antibody binds human factor XII in solid-phase ELISA and immunoblotting techniques.

Approximate working concentration is 1-5 ug/mL. Researcher should titer antibody in specific assay.

Reconstitute with deionized water, aliquot and store at -20C for prolonged periods. Avoid freeze-thaw cycles.

Thank You,
Quality Assurance

Sample