



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Lot:** SAMPLE SAMPLE

**Description:** Monoclonal antibody to Human Protein C

**Concentration:** 1 mg/ml (before lyophilization)

**Aliquots:** 1 x 0.50 mg

**Volume/aliquot:** 0.500 ml

**Activity:** N/A

**Total protein:** 0.50 mg

**Buffer:** 10 mM Sodium Phosphate/ 0.15 M  
NaCl/ 0.01% Sodium Azide/ pH 7.4

**Storage:** -20 °C

**Expiration Date:** 01/2010

**PO:** N/A

Immunogen: Human Protein C

Clone ID#: 067

Purified from cell supernatant by Protein G affinity chromatography.

Binds human Protein C in solid-phase ELISA. Blots non-reduced and heavy chain of reduced human Protein C.

We suggest hydrating the antibody to original volume with water.

Thank You,  
**Quality Assurance**

Sample