



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Plasminogen

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 10/2010

PO: N/A

Clone ID#: 039

Antibody binds human plasminogen in solid-phase and solution-phase ELISA, and human plasminogen in immunoblotting techniques. Antibody does not bind plasminogen kringle fragments or recombinant angiostatin.

Reconstitute to original volume with water. Avoid repeated freeze-thaw cycles.

Thank You,
Quality Assurance

Sample