



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to human Plasminogen and Angiostatin Kringle 1-4

Concentration: 1 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 05/2009

PO: N/A

Antibody binds recombinant human angiostatin and human plasminogen, specifically kringles 1-4, in solid-phase and solution-phase ELISA. Antibody binds only slightly to isolated plasminogen kringle fragment 1-3 and isolated fragment 4. It binds human plasminogen and recombinant angiostatin in immunoblotting techniques.

Recommended use: 1-5 ug/mL. Researcher should titer antibody in specific assay.

Store lyophilized antibody at -20C. Reconstitute with good laboratory water and centrifuge to remove insoluble material. Store reconstituted antibody in aliquots at -20C. Avoid repeated freeze-thaw cycles.

Thank You,
Quality Assurance

Sample