



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Monoclonal antibody to Human Vitronectin

Concentration: 1.0 mg/ml (before lyophilization)

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.500 ml

Activity: N/A

Total protein: 0.50 mg

Buffer: 10 mM Sodium Phosphate/ 0.15 M
NaCl/ 0.01% Sodium Azide/ pH 7.4

Storage: -20 °C

Expiration Date: 06/2009

PO: N/A

Antibody purified by Protein G affinity chromatography from cell supernatant. Purity determined by SDS-PAGE to be >98%.
Antibody binds human vitronectin in solid phase ELISA and immunoblotting procedures. Binds plasma vitronectin when coupled to sepharose. Recommended use: 1-5 ug/mL. Reconstitute with water to original concentration.

Thank You,
Quality Assurance

Sample