



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: PIIa

Lot: SAMPLE SAMPLE

Description: Porcine Thrombin

Concentration: 1.80 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.056 ml

Activity: 3047.00 NIH u/mg

Total protein: 0.10 mg

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/ 0.1% PEG-8000/ pH 6.5

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Porcine Thrombin was activated from homogeneous porcine prothrombin by activation with Bovine Factor Xa, Factor Va, and phospholipid. These activating enzymes were removed after activation.

Thank You,
Quality Assurance

Sample