



**1801 Commerce Drive**  
**South Bend, IN 46628**  
**Phone:(574) 288-2268**  
**Fax:(574) 288-2272**  
**Email: info@enzymeresearch.com**

## Certificate of Analysis

**Catalog #:** RBX

**Lot:** SAMPLE SAMPLE

**Description:** Rabbit Factor X

**Concentration:** 1.93 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.052 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ 1mM Benzamidine/ pH 7.4

**Storage:** < -60° C

**Expiration Date:** SAMPLE SAMPLE

**PO:** N/A

**Date of Manufacture:** SAMPLE SAMPLE

Rabbit Factor X is purified via salt precipitations along with ion and affinity chromatography. It is judged >94% pure on SDS PAGE.

Thank You,  
**Quality Assurance**

Sample