



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com

Certificate of Analysis

Catalog #: RBXa

Lot: SAMPLE SAMPLE

Description: Rabbit Factor Xa

Concentration: 1.71 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.058 ml

Activity: 171.57 IU/mg

Total protein: 0.10 mg

Buffer: 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

Rabbit Factor Xa is prepared from homogeneous Rabbit Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation is observed via SDS-PAGE.

Thank You,
Quality Assurance

Sample