



**1801 Commerce Drive**  
**South Bend, IN 46628**  
**Phone:(574) 288-2268**  
**Fax:(574) 288-2272**  
**Email: info@enzymeresearch.com**

## Certificate of Analysis

**Catalog #:** RFXa

**Lot:** SAMPLE SAMPLE

**Description:** Rat Factor Xa

**Concentration:** 1.52 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.066 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 20 mM Tris-HCl/ 0.7 M NaCl/ pH 7.4

**Storage:** < -60° C

**Expiration Date:** SAMPLE SAMPLE

**PO:** N/A

**Date of Manufacture:** SAMPLE SAMPLE

Rat Factor Xa is prepared from homogenous Rat Factor X by activation with Russells' Viper Venom. This RVV-X is removed after activation. Complete activation is observed on SDS-PAGE.

Thank You,  
**Quality Assurance**

Sample