



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: R11a

Lot: SAMPLE SAMPLE

Description: Rat Thrombin

Concentration: 1.26 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.079 ml

Activity: 3033.00 NIH u/mg

Total protein: 0.10 mg

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/ 0.1% PEG-8000/ pH 6.5

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The rat thrombin was activated from plasma purified rat prothrombin by activation with venom. This venom is removed after activation. Material is >90% pure by SDS-PAGE.

Thank You,
Quality Assurance

Sample