



**1801 Commerce Drive  
South Bend, IN 46628  
Phone:(574) 288-2268  
Fax:(574) 288-2272  
Email:**

## Certificate of Analysis

**Catalog #:** APCi

**Lot:** Sample

**Description:** Activated Protein C Inactivated

**Concentration:** 0.79 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.127 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 20 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

**Storage:** < -60 C

**Expiration Date:** Sample

**PO:** Sample

**Manufacture Date:**

Our Activated Protein C is inactivated using PPACK, Dihydrochloride and is removed via gel filtration. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**