



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: FN

Lot: SAMPLE SAMPLE

Description: Fibronectin

Concentration: 1.44 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.694 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 100mM Hepes/ 150mM NaCl/ 5% Sucrose/ pH 7.4

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

Fibronectin is purified from fresh frozen plasma via gelatin sepharose. It is a glycoprotein, composed of two subunits of approximately 220 kD. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV DNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample