



**1801 Commerce Drive**  
**South Bend, IN 46628**  
**Phone:(574) 288-2268**  
**Fax:(574) 288-2272**  
**Email: info@enzymeresearch.com**

## Certificate of Analysis

**Catalog #:** HAPO-H

**Lot:** SAMPLE SAMPLE

**Description:** Human Apolipoprotein-H (B2 Glycoprotein-1)

**Concentration:** 1.60 mg/ml

**Aliquots:** 1 x 0.10 mg

**Volume/aliquot:** 0.062 ml

**Activity:** N/A

**Total protein:** 0.10 mg

**Buffer:** 20 mM EPPs/ 0.15 M NaCl/ pH 8.0

**Storage:** < -60° C

**Expiration Date:** SAMPLE SAMPLE

**PO:** N/A

**Date of Manufacture:** SAMPLE SAMPLE

The above protein was purified from human plasma via ion and affinity column chromatography. It is >95% pure on SDS-PAGE.

The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**

Sample