



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com

Certificate of Analysis

Catalog #: HFIXa 1080

Lot: SAMPLE SAMPLE

Description: Human Factor IXa beta

Concentration: 1.12 mg/ml

Aliquots: 1 x 0.50 mg

Volume/aliquot: 0.446 ml

Activity: 6.85 U/mg

Total protein: 0.50 mg

Buffer: 50mM MES / 0.6M NaCl / pH 6.0 -NaOH

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Human Factor IXab was activated from homogeneous Human Factor IX by activation with Human Factor XIa. The XIa is removed after activation. Complete activation is observed on SDS-PAGE. The above protein was purified from Human plasma that was tested and found negative by FDA accepted methods for Anti-HIV1/2, Anti-HTLV I & II, HBsAg, Anti-HCV, Syphilis, HBC Ab, HIV-1 p24 Ag or HIV-1 RNA, HCV RNA and HBV RNA. Donors are screened for CJD (Creutzfeldt-Jakob Disease).

Thank You,
Quality Assurance

Sample