



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HFX-GD

Lot: Sample

Description: Human Factor X-GD

Concentration: 1.71 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.585 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 100 mM Tris-HCl/ 50 mM NaCl/ pH 7.5

Storage: < -60 C

Expiration Date: Sample

PO: Sample

Date of Manufacture:

Factor X- gla domainless was prepared from purified Human Factor X using alpha-chymotrypsin. The alpha-chymotrypsin and all other products were removed using an anion exchange column after digestion. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance