



**1801 Commerce Drive**  
**South Bend, IN 46628**  
**Phone:(574) 288-2268**  
**Fax:(574) 288-2272**  
**Email:**

## Certificate of Analysis

**Catalog #:** HFXabGD

**Lot:** Sample

**Description:** Human Factor Xab gla domainless

**Concentration:** 2.05 mg/ml

**Aliquots:** 1 x 1.00 mg

**Volume/aliquot:** 0.488 ml

**Activity:** 191.00 IU/mg

**Total protein:** 1.00 mg

**Buffer:** 20mM Bis-Tris/ 0.1M NaCl / pH 6.0

**Storage:** < -60 C

**Expiration Date:** Sample

**PO:** Sample

**Manufacture Date:**

The Human Factor Xa beta gla domainless was prepared from HFX-GD using purified RVV-X. This RVV-X was removed after activation by gel filtration. Complete activation was observed on SDS PAGE. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,  
**Quality Assurance**