



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email:**

Certificate of Analysis

Catalog #: HK 1300

Lot: Sample

Description: Single Chain HMW Kininogen

Concentration: 1.91 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.524 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 4 mM sodium acetate-HCl/ 0.15 M NaCl/ pH 5.3

Storage: < -60 C

Expiration Date: Sample

PO:

Date of Manufacture: Sample

The Single Chain High MW Kininogen is greater than 95% single chain. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance