



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com

Certificate of Analysis

Catalog #: HPG 2001

Lot: SAMPLE SAMPLE

Description: Human Glu-Plasminogen

Concentration: 3.12 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.321 ml

Activity: 20.50 CU/mg

Total protein: 1.00 mg

Buffer: 50 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Human glu-plasminogen is homogeneous on SDS-PAGE gels and shows no reduction upon incubation with 2-mercaptoethanol. No Plasmin activity was detected using the chromogenic substrate S-2251. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance

Sample