



**1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com**

Certificate of Analysis

Catalog #: LPG 2002

Lot: SAMPLE SAMPLE

Description: Human Lys-Plasminogen

Concentration: 0.97 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 1.031 ml

Activity: N/A

Total protein: 1.00 mg

Buffer: 50 mM Tris-HCl/ 0.1 M NaCl/ pH 7.4

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Human Lys-Plasminogen was prepared from homogeneous glu-plasminogen by batching with a plasmin-sepharose column. This resin was removed after activation. Complete activation was observed by SDS-PAGE. No plasmin activity was detected using the chromogenic assay S-2251. The above protein was purified from Human plasma units that test negative for antibodies to HBc, HCV, HIV 1/2, HTLV-I/II, and non reactive for HBsAg, HIV-1 RNA, HCV RNA, HBV DNA and Syphilis by FDA approved tests. Donors are screened for CJD.

Thank You,
Quality Assurance

Sample