



1801 Commerce Drive
South Bend, IN 46628
Phone:(574) 288-2268
Fax:(574) 288-2272
Email: info@enzymeresearch.com

Certificate of Analysis

Catalog #: RFib

Lot: SAMPLE SAMPLE

Description: Rat Fibrinogen

Concentration: 23.49 mg/ml

Aliquots: 1 x 1.00 mg

Volume/aliquot: 0.043 ml

Activity: 94.40 % Clottable

Total protein: 1.00 mg

Buffer: 20 mM sodium citrate-HCl/ pH 7.4

Storage: < -60° C

Expiration Date: SAMPLE SAMPLE

PO: N/A

Date of Manufacture: SAMPLE SAMPLE

The Rat Fibrinogen was purified from fresh frozen rat plasma using a combination of salt precipitations and column chromatography. The purified protein is homogeneous on SDS-PAGE and shows greater than 90% clottability in a thrombin based assay.

Fibrinogen MUST be thawed in a 37C water bath. Once thawed, do not chill. Remaining protein may be aliquoted into single use size and frozen at <-60C.

Thank You,
Quality Assurance

Sample